




1. Insert band through outer slot



2. Return band through middle slot



3. Insert band through inner slot



**4. Return band through outer slot.
5. Repeat on other side**

Clean gently using a soft cloth with soap & water

IMPORTANT DISCLAIMER

This product has not been tested or evaluated for the prevention of any specific diseases or infections, the filtering of specific amounts of viruses or bacteria, or the elimination or reduction of the amount of viruses or bacteria. Individual users of this shield are completely responsible for their own decisions regarding appropriate personal protective equipment for themselves. The vendor and manufacturer of this product provides no covenant, representations or warranties that the product can protect against any virus or bacteria or any illness and use of the product is at your own risk. Any claim and damages related thereto against the vendor or manufacturer shall in any case be limited to the cost of the product. Always insure that the product is cleaned and disinfected properly after each use. Use of this product shall be deemed to be acknowledgement and acceptance of these terms and conditions and limits of liability.

In an effort to mitigate the spread of coronavirus disease please remember to practice social distancing, wash your hands often with soap, use alcohol-based hand sanitizer, avoid touching your eyes, nose and mouth, cough and sneeze into your elbow, use cleaning products to disinfect objects and surfaces you commonly touch and stay at home if you are feeling sick, as recommended by the Public Health Agency of Canada.

Scan for online instruction videos



1 800 410-0264 - displaypeople.com